



ORANKA system supplements.
Clever solutions for your individual needs.



ORANKA offers additional clever system solutions suitable for its award-winning fruit drink concentrates and practical drinks dispensers. For example, the ORANKA 2.5-litre carafe system comes in two variants. The plastic carafe fits in any fridge door and offers you outstanding handling and practical markings for precise dispensing. The extremely high-quality glass carafe with cooling rod can also be equipped with a second cooling rod as an option

The very handy ORANKA 1.1-litre glass carafe is, of course, also available in dishwasher-safe quality. The easy handling and space-saving dimensions also allow practical temporary storage in the fridge. For easy dispensing, we also offer dishwasher-safe glasses and cups.

ORANKA system enhancements.

Clever solutions for your individual needs.



Glasses and cups

- ① 6 Drinking glasses, Standard, 0.2 litres (item no. 9002000)
- ② 6 Drinking glasses, Exclusive, 0.2 litres (item no. 9002003)
- ③ Drinking cup, 0.33 litres with lid (item no. 9002010)
- ④ Drinking cup, 0.33 litres without lid (item no. 9002011)
- ⑤ Lid for drinking cups, 0.33 litres without lid (item no. 9002012)

Karaffen und Zubehör

- ⑥ Serving tray ø 37 cm, ORANKA Catering-Pro (item no. 9005000)
With anti-slip coating (Hand wash recommended)
- ⑦ Glass carafe, 2.5 litres with cooling rod (item no. 9004002)
- ⑧ 6 glass carafes, 2.5 litres with cooling rod in box (item no. 9004000)
- ⑧ Cooling rod for glass carafe, 2.5 litres (item no. 9004001)
- ⑨ Glass carafe, 1.1 litres with lid (item no. 9024000)
- ⑩ Plastic carafe, 2.5 litres with lid (item no. 9034000)

Measuring cups and dispensing containers

- ⑪ Measuring cup, 5 litres (item no. 9305000)
- ⑫ Dispensing container, 10 litres (item no. 9310000)



Free delivery when you order concentrate at the same time!

All prices subject to change, plus 19 % VAT. Delivery: in accordance with our Terms & Conditions. Payment terms: by individual agreement.



Practical dispensing aid

With the handy printed markings, ORANKA concentrates can be dispensed precisely and prepared directly in the carafe with no complicated calculations.

